


ClearViewTM

Anti-Beschlag Optik Wärmer



Heizt die Optik auf Körpertemperatur vor

Kein Strom oder heißes Wasser

Aktiviert durch einfaches Anklicken



supporting surgical outcomes

ClearView™


Anti-Beschlag Optik Wärmer



KEBOMED

Anwendung

1. OPTIK IN CLEARVIEW LEGEN



2. SCHEIBE BIEGEN (2-3 X)




3. ERWÄRMT FÜR DEN GEBRAUCH

Achtung: Der Optikwärmer sollte 15 Minuten vor Gebrauch der Optik aktiviert werden.

ClearView™ bleibt bis zu 2 Stunden warm. Durch Abdecken mit einem Tuch kann diese Zeit verlängert werden.

References

Clearview scopewarmer generates a significant higher difference during the complete process of pre-warming the endoscope. Clearview exchanges 75% more energy at the start compared to a warm water bath. Over time the diffusion capacity of a warm bath reduces to only 12% of the heat flux from the Clearview scopewarmer*. This significant difference ensures that enough energy is diffused during the process to ensure the whole endoscope, including the optic materials inside, are pre-warmed to body temperature.



* Data on file

Für Deutschland:
KEBOMED GAS AG
PHONE +41 41 450 3410
FAX +41 41 450 3411
EMAIL: INFO@KEBOMED.CH

KEBOMED